



WWW.DENTALLEARNING.NET

# DENTAL LEARNING

A PEER-REVIEWED PUBLICATION

*Knowledge for Clinical Practice*



## Single-Patient-Use VS Multiple-Patient-Use Dental Burs

*By Keith Cyrus Keller, DMD*

INSIDE  
Earn 2  
**CE**  
Credits

Written for  
dentists, hygienists,  
and assistants

**ADA CERP**® | Continuing Education  
Recognition Program

Integrated Media Solutions Inc./DentalLearning.net is an ADA CERP Recognized Provider. ADA CERP is a service of the American Dental Association to assist dental professionals in identifying quality providers of continuing dental education. ADA CERP does not approve or endorse individual courses or instructors, nor does it imply acceptance of credit hours by boards of dentistry. Concerns or complaints about a CE provider may be directed to the provider or to ADA CERP at [www.ada.org/cerp](http://www.ada.org/cerp). Integrated Media Solutions Inc./Dental Learning.net designates this activity for 2 continuing education credits.



Approved PACE Program Provider  
FAGD/MAGD Credit Approval  
does not imply acceptance by a  
state or provincial board of  
dentistry or AGD endorsement.  
2/1/2020 – 1/31/2024  
Provider ID: # 346890  
AGD Subject Codes: 148, 250



Dental Learning, LLC is a Dental Board of California CE Provider. The California Provider # is RP5062. All of the information contained on this certificate is truthful and accurate. Completion of this course does not constitute authorization for the attendee to perform any services that he or she is not legally authorized to perform based on his or her license or permit type. This course meets the Dental Board of California's requirements for 2 units of continuing education. CA course code is 02-5062-20013



## ABSTRACT

Dental armamentarium has evolved since the first tool-assisted manipulation of teeth thousands of years ago. In the modern era, discussion continues as to which tools will provide the best outcomes for patients and dental providers. One significant debate, especially in light of the heightened awareness of infection control, pits dental burs that are single-patient-use vs multiple-patient-use. Important considerations include overhead costs, cutting efficiency, and infection control. Manufacturing of dental burs and other armamentarium has an interesting history, and in this course you will learn where the technology and science stand at this moment in time.

## EDUCATIONAL OBJECTIVES

After completing this course, the reader should be able to:

- Discuss the history of dental burs and handpieces
- List the pros and cons of single-patient-use and multiple-patient-use armamentarium
- Analyze dental bur cost, efficiency, and infection control.

## ABOUT THE AUTHOR



By Keith Cyrus Keller, DMD

Keith Cyrus Keller, DMD, graduated magna cum laude from Tufts University School of Dental Medicine. After dental school, he completed a general practice residency at Palmetto Health Richland Hospital in Columbia, SC, and three years of training in orthodontics and dentofacial orthopedics at St. Barnabas Hospital in Bronx, NY. Currently, Dr. Keller is an attending at Bronx Care hospital, where he lectures and trains residents studying orthodontics. He has a private practice in Brooklyn, NY.

## Introduction

Removing tooth structure is an essential procedure in dentistry. Archeological evidence (Figure 1) from 13,000 years ago in Tuscany, Italy, has shown “tool-assisted manipulation to remove necrotic or infected pulp in vivo and the subsequent use of a composite, organic filling.”<sup>1</sup> Other findings include 11 drilled molar crowns from nine adults discovered in a graveyard in Pakistan that dates from 7,500 to 9,000 years ago.<sup>2</sup>

For thousands of years, hand-held drills could produce roughly 15 rotations per minute (rpm). In

SPONSOR/PROVIDER: This is a Dental Learning, LLC continuing education activity. COMMERCIAL SUPPORT This educational activity is made possible through an unrestricted educational grant from Microcopy. STATEMENTS: Dental Learning, LLC is an ADA CERP recognized provider. ADA CERP is a service of the American Dental Association to assist dental professionals in identifying quality providers of continuing dental education. ADA CERP does not approve or endorse individual courses or instructors, nor does it imply acceptance of credit hours by boards of dentistry. Dental Learning, LLC designates this activity for 2 CE credits. Dental Learning, LLC is also designated as an Approved PACE Program Provider by the Academy of General Dentistry. The formal continuing education programs of this program provider are accepted by AGD for Fellowship, Mastership, and membership maintenance credit. Approval does not imply acceptance by a state or provincial board of dentistry or AGD endorsement. The current term of approval extends from 2/1/2020 – 1/31/2024. Provider ID: # 346890. EDUCATIONAL METHODS: This course is a self-instructional journal and web activity. Information shared in this course is based on current information and evidence. REGISTRATION: The cost of this CE course is \$29.00 for 2 CE credits. PUBLICATION DATE: April 2020. EXPIRATION DATE: March 2023. REQUIREMENTS FOR SUCCESSFUL COMPLETION: To obtain 2 CE credits for this educational activity, participants must pay the required fee, review the material, complete the course evaluation and obtain a score of at least 70%. AUTHENTICITY STATEMENT: The images in this course have not been altered. SCIENTIFIC INTEGRITY STATEMENT: Information shared in this continuing education activity is developed from clinical research and represents the most current information available from evidence-based dentistry. KNOWN BENEFITS AND LIMITATIONS: Information in this continuing education activity is derived from data and information obtained from the reference section. EDUCATIONAL DISCLAIMER: Completing a single continuing education course does not provide enough information to result in the participant being an expert in the field related to the course topic. It is a combination of many educational courses and clinical experience that allows the participant to develop skills and expertise. PROVIDER DISCLOSURE: Dental Learning does not have a leadership position or a commercial interest in any products that are mentioned in this article. No manufacturer or third party has had any input into the development of course content. CE PLANNER DISCLOSURE: The planner of this course, Mary Benedon, does not have a leadership or commercial interest in any product or services discussed in this educational activity. Any questions or comments can be sent to Support@dentallearning.net. TARGET AUDIENCE: This course was written for dentists, dental hygienists, and assistants, from novice to skilled. CANCELLATION/REFUND POLICY: Any participant who is not 100% satisfied with this course can request a full refund by contacting Dental Learning, LLC in writing or by calling 1-888-724-5230. Please direct all questions pertaining to Dental Learning, LLC or the administration of this course to Support@dentallearning.net. Go Green, Go Online to www.dentallearning.net to take this course. © 2020

Copyright 2020 by Dental Learning, LLC. No part of this publication may be reproduced or transmitted in any form without prewritten permission from the publisher.



500 Craig Road, First Floor, Manalapan, NJ 07726

President  
ALDO EAGLE

CE Project Manager  
MARY BENEDON

Creative Director  
MICHAEL HUBERT

Art Director  
JOE CAPUTO



1864, British dentist George Fellows Harrington invented a mechanical, clockwork drill called Erado (from the Latin for “I scrape out”). It was the first continuously rotating drill that lasted 2 minutes (Figure 2). James Morrison devised a pedal-powered drill (2,000 rpm) in 1871; it operated on the same principle as the treadle sewing machine (Figure 3). George Green revolutionized dentistry in 1875 with a patent for the first electric drill (Figure 4). By 1914, electric drills could reach speeds of up to 3,000 rpm. A second wave of rapid development occurred in the 1950s and 1960s, including the development of the air-turbine drill. Current iterations can operate at up to 800,000 rpm; however, 400,000 rpm (high speed) and 40,000 rpm (low speed) handpieces are most common.


In addition to drills, dentists use metal, hand-held instruments in a variety of shapes that act as miniature chisels and excavators to remove caries and alter the architecture of a tooth in preparation for a restorative material. Although new techniques and technologies—such as lasers, air abrasion, and chemical dissolution—are advocated for removing dental hard tissue during tooth preparation, high- and low-speed rotary instruments remain the workhorses of modern dentistry and are used for caries removal, gross tooth reduction, margin refinement, enameloplasty, restoration removal, surgical extractions, crown lengthening, orthodontic adhesive removal, implant placement, and acrylic and metal appliance adjustment.

## Dental Burs

The wide variety of available dental burs helps dentists accomplish many daily tasks (Figure 5). The International Organization for Standardization introduced a 15-digit code numbering system in 1979. The most up-to-date standard (ISO 6360-1:2004) helps to identify general and specific characteristics:

- The first group of three digits identifies materials used for the working part of instruments.
- The second group of three digits identifies instrument shanks and handles and overall instrument length.
- The third group of three digits identifies instrument shapes.
- The fourth group of three digits identifies specific characteristics for groups of instruments.
- The fifth group of three digits identifies the nominal diameter of the working part of the instruments (nominal size).<sup>3</sup>

The working part of a modern dental bur can be made of stone, stainless steel, stainless steel coated with tungsten carbide, pure tungsten carbide, diamond particles, ceramic, or plastic.



**Figure 1**—The upper right and left first incisors (RI1, LI1) of Riparo Fredian. (A) The RI1 in occlusal view. (B) Volume rendering of RI1 viewed from the lingual side. (C) Digital reconstruction of the RI1 with transparent dentine to show the residue; mesio-lingual view. (D) The LI1 in occlusal view; note the black patina within the cavity. (E) Volume rendering of LI1 viewed from the lingual side. (F) Digital reconstruction of LI1 with transparent dentine to show the residue: distal view. Scale bar: a and d 2 mm; b, c, e, and f 1 cm.





## Tungsten Carbide


In chemistry, “carbide” describes a compound composed of carbon and a metal. Tungsten carbide contains equal parts tungsten (Figure 6) and carbon atoms. In its most basic form, tungsten carbide is a fine gray powder, but it can be pressed and formed into shapes via sintering. Compared to steel, tungsten carbide is approximately twice as stiff and double the density.<sup>4</sup>

In one study, five experienced operators were given a set of No. 330 tungsten carbide burs that had been separated into five blind categories: those not autoclaved at all and those that had been autoclaved one, two, five, or 10 times. Each operator was given a freshly extracted tooth and told


to judge the cutting effectiveness of each bur by making an initial plunge into the tooth and then completing an occlusal cavity preparation in molar models made from composite. The results indicate that autoclaving tungsten carbide burs doesn’t affect the subjective perception of effectiveness.<sup>5</sup>

## Diamond

The first diamond burs were invented in 1897 by Willman and Schroeder from the University of Berlin, Germany.<sup>6</sup> These early burs were made by hammering diamond powder into the surface of soft copper or iron blanks. Silicone carbide burs had been used before this time, but they were found to wear when used on enamel. The modern diamond bur was created in 1932 by German industrialist W.H. Drendel, who developed a process for bonding diamond points to stainless steel.



**Figure 2**—George Fellows Harrington's first clockwork dental drill in a silver-plated veneered box. Invented in 1864, it had interchangeable drill heads and was operated by a spring-wound mechanism wound with a key.



**Figure 3**—Dental foot engine designed by James B. Morrison and manufactured by S.S. White Dental Manufacturing Company.



The most common bonding method today is electrolytic co-deposition of natural or synthetic diamond particles within a nickel or duplex nickel–chromium matrix onto a stainless steel shank.<sup>7</sup> The diamond particles are embedded into the metal matrix, like chopped nuts embedded in caramel that coats a candied apple (Figure 7). Other techniques to attach diamond particles include microbrazing, chemical vapor deposition, sintering, and adhesives.<sup>8</sup> A survey of 131 dentists showed a preference

for diamond burs (75%) over tungsten carbide burs (15%) and that most (85%) used their burs until they were worn out. Shockingly, clinicians usually sterilized burs either once daily (35%) or for every patient (35%).<sup>9</sup>

## Cutting Efficiency

In a study by Pilcher and colleagues, a scanning electron microscope was used to look at six different single-patient-use and two different multiple-patient-use diamond burs after making 20 cuts into a block of glass ceramic. Gradual reduction in cutting rates was observed and all eight burs exhibited diamond particle wear, as well




Figure 4—1875 patent of first electric drill



Figure 5—Burs with a variety of shapes, sizes, and materials



Figure 6—Tungsten rods with evaporated crystals, partially oxidized with colorful tarnish. Purity 99.98%, as well as a high pure (99.999% = 5N) 1 cm<sup>3</sup> tungsten cube for comparison.



as loss of diamond particles and binder. Five of six single-patient-use burs had mean volumetric cutting rates similar to the multiple-patient-use diamonds.<sup>10</sup>

Another study looked at 16 different burs (some single-use and some multipatient-use) from a variety of manufacturers. “The results showed no differences in cutting efficiency between medium-grit conventional and disposable diamond burs.”<sup>11</sup>

How often can a multiple-patient-use diamond be used? Researchers concluded that “diamond burs wear after multiple use and that they should be changed after

five teeth preparations at most. A diamond bur should not be used for teeth preparation after try-in procedures of metal or zirconia substructures.”<sup>12</sup>

What is happening at a microscopic level? Several theories have been given for the decreased efficiency of diamond burs, including diamond wear-out, diamond pull-out, clogging by debris, and degradation of the diamond binder material. In a study of 15 new diamond burs, a scanning electron microscope was used to view the same area of the bur before and after use on human teeth. The results showed that diamond pull-out did not contribute more than 10% to the total wear. Clogging was responsible for the wear of two out of 15 burs. Particle wear-out, with blunting of the sharp edges and intensive diamond mass loss, was found to be the dominant wear mechanism. Damage to the binder layer was not detected.<sup>13</sup>

### Overhead Cost

Modern production techniques have resulted in a decrease in the price-point comparisons of tungsten carbide and diamond burs. Manufacturers offer single-patient-use burs that are less expensive than multiple-patient-use burs, but are they cost effective? Three aspects should be considered: initial bur cost, reusable bur longevity, and overhead costs (including the non-income-generating labor involved in cleaning, sorting, sterilizing, storing, and retrieving burs).

In a study comparing five different single-use and five different multiple-patient-use diamond burs, single-use burs were found to be more cost effective. When analyzing the longevity of the burs to perform full crown preps on lower first molars, all 10 brands cut the first crown within 2 minutes. All 10 brands cut more slowly on each subsequent preparation, and cut times doubled for all brands after about four preps. Overall, the author concluded that single-patient-use diamonds used for one crown prep had the best combination of speed and low




Figure 7—SEM micrographs showing (a) the morphology of a step diamond bur 12S, (b) diamond abrasives on the bur with approximately 64-µm grit size (Alao et al., 2017). Reprinted with permission from Elsevier.



Figure 8—Microcopy single-patient-use diamond burs.

cost.<sup>14</sup> As noted previously, diamond burs should not be used for more than five tooth preparations. Following this recommendation, the price comparison favors single-patient-use burs (considering initial cost only).

It is time-consuming for staff to properly follow protocol when attempting to reuse burs. In addition to the time required and the questionable efficacy of reprocessing, personnel also must be aware of the potential for accidental sharps exposure, especially when hand scrubbing is used for presterilization cleaning. In a survey of 164 dentists, almost a third (51) reported having had a sharp injury within their career, needlestick (32%) and dental bur (26%) accidents were the most common.<sup>15</sup>

### Infection Control

Three sterilization methods to remove bacterial spores, bacteria, viruses, molds, and fungi currently are used in dentistry: steam under pressure in a steam autoclave, dry heat in a sterilizing oven, and chemical vapor.<sup>16</sup> Each method is technique sensitive and has pros and cons. Are routine dental sterilization procedures effective? In one study, previously used burs and endodontic files were gathered from different dental offices after they had been packaged and sterilized for reuse. Sterilization procedures were carried out according to the protocols used by each office and, surprisingly, each technique was less than 100% effective.<sup>17</sup> The problem is the method used to remove gross debris before sterilization. However, if proper protocols are followed, most organic debris can be physically removed, making it theoretically possible to sterilize instruments. But which protocol is most effective?

Experiments conducted using *Streptococcus sanguinis*, a pathogen identified as the causative organism in bacterial endocarditis and meningitis,<sup>18</sup> have shown that autoclaving alone was not effective in decontaminating steel or diamond burs. The authors then tested three presterilization cleaning techniques: manual scrubbing,

soaking in an enzymatic bath, and washer-disinfectors. Manual scrubbing, which was shown to be operator sensitive and inconsistent, is not recommended due to serious risk of cross-infection from puncture wounds. Enzymatic agents were very effective, but their use also is operator-dependent. The researchers concluded: “The washer-disinfectors investigated in this study rendered all test bacteria non-viable. Burs contaminated in practice may harbor spore-forming bacteria and blood-borne viruses that may not be killed by a washer-disinfectant and should be autoclaved after precleaning.” Due to the reproducibility, ease of operation, and minimal contact with contaminated instruments, washer-disinfectors were the most effective method of presterilization cleaning. “It is suggested that enzymatic agents should be confined for use immediately postoperatively, for immersion of instruments prior to washer-disinfectors.”<sup>19</sup>

Of particular concern is the iatrogenic transmission of rare, fatal degenerative disorders of the nervous system, such as Creutzfeldt–Jakob disease, Gerstmann–Straussler–Scheinker disease, Kuru, and fatal familial insomnia. These diseases are caused by prions, infectious agents smaller than viruses that contain no DNA or RNA. Stanley B. Prusiner, MD, distinguished these infectious particles from viruses or viroids and designated them as prion proteins (PrP); in 1997, he was awarded the Nobel Prize in Physiology or Medicine. The incubation period ranges from two to 35 years, which means that people may transmit the diseases while asymptomatic for decades.<sup>20</sup> Prions are far more resistant to physical and chemical inactivation than even bacterial spores. The only sterilization methods that appear to be completely effective are sodium hypochlorite solutions or hot solutions of sodium hydroxide.<sup>21</sup> World Health Organization guidelines suggest the use of immersion in sodium hypochlorite (20,000 ppm available chlorine) for 1 hour, boiling in 1 mole/liter of sodium hydroxide





for 1 hour, or vacuum autoclaving at 121°C for 30 to 90 minutes in the presence of 2 moles/liter of sodium hydroxide.<sup>21</sup> “Such procedures are inappropriate for surgical or dental instruments due to degradation and tarnishing of the metal surface, hazards to the user and issues with disposal.”<sup>22</sup> As of 2007, the United Kingdom has banned the reuse of endodontic files.

It has become common practice to minimize infectious particle transmission by using single-use armamentarium: masks, gloves, anesthetic carpules, syringe needles, scalpel blades, cotton rolls, gauze, prophylaxis cups, matrix bands, dental dams, saliva ejectors, impression trays, air–water syringe tips, bibs, and headrest and seat covers. The workflow required to reuse burs places an additional burden on staff to control the spread of infection from staff to patient, patient to staff, and patient to patient. Additionally, patients have peace of mind when they experience a well-run office that demonstrates infection-control techniques, including opening sterile packaging at chairside. For example, Microcopy manufactures a variety of diamond (Figure 8) and tungsten carbide burs designed for single-patient-use. Following the ISO 11137, which governs the sterilization of healthcare products, each bur is identified with “STERILE R” to indicate it is sterilized by gamma radiation.<sup>23,24</sup>

## Centers for Disease Control and Prevention Recommendations

“Sometimes disposable patient-care items have reusable heat-tolerant alternatives. This is often true for prophylaxis angles, high-volume evacuator tips, impression trays, dental burs, and air and water syringe tips. In many instances, the reusable version of these items is hard to clean adequately, and it may be safer, easier, and more cost-effective to use the disposable version. When determining the cost-effectiveness of disposable and reusable items, dental healthcare personnel should consider not only the cost of

the disposable item, but also the cost, time, and materials involved with cleaning and reprocessing the reusable item.

“Some devices—such as burs, endodontic files, and broaches—may be practical to consider single-use because the way they are constructed makes them hard to clean. In addition, cleaning and heat sterilization can lead to deterioration on the cutting surfaces and raise the potential for breakage during patient treatment. The Food and Drug Administration (FDA) considers all diamond-coated burs and scaler tips single-use unless their manufacturers have submitted a 510(k) for reprocessing...Dental healthcare personnel (DHCP) should always refer to manufacturer instructions to determine if a device is single-use. If a device does not have validated reprocessing (i.e., cleaning and disinfection or sterilization) instructions, it is considered single-use (i.e., disposable).”<sup>25</sup>

“During aerosol-generating procedures...in patients who are not suspected of being infected with an agent for which respiratory protection is otherwise recommended (e.g., M. tuberculosis, SARS, or hemorrhagic fever viruses), wear one of the following: a face shield that fully covers the front and sides of the face, a mask with attached shield, or a mask and goggles (in addition to gloves and gown).”<sup>26</sup>

## Conclusion

Cost and efficiency are significant factors when operating a thriving practice, but infection control is paramount. Regulatory frameworks exist mainly in an effort to control infection with certain items mandated as single-patient-use. Infection control has become a focal point in the modern era, and the decision to move away from multiple-patient-use to single-patient-use armamentarium is ultimately up to the dentist. There has been an evolution in dental technology and scientific understanding of disease transmission, with concerns that are not unfounded. The AIDS pandemic revolutionized the use of personal protective equipment in the 1980s, and with the recent





outbreak of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and the associated coronavirus disease of 2019 (COVID-19), the dental world is poised to further embrace single-patient-use armamentarium and respirators during aerosol-producing procedures.

## References

1. Oxilia G, Fiorillo F, Boschin F, et al. The dawn of dentistry in the late upper Paleolithic: An early case of pathological intervention at Riparo Fredian. *Am J Phys Anthropol*. 2017;163(3):446-61.
2. Coppa A, Bondioli L, Cucina A, et al. Palaeontology: Early Neolithic tradition of dentistry. *Nature*. 2006;440(7085):755-6.
3. International Organization for Standardization. ISO 6360-1:2004. Dentistry—Number coding system for rotary instruments—Part 1: General characteristics. April 2004. [iso.org/standard/37301.html](http://iso.org/standard/37301.html)
4. Wikipedia. Tungsten carbide. [en.wikipedia.org/wiki/Tungsten\\_carbide](http://en.wikipedia.org/wiki/Tungsten_carbide)
5. Spranley TJ, Cheramie TJ, Ireland EJ, et al. Cutting effectiveness of carbide burs following multiple steam autoclaving cycles. *Gen Dent*. 2011;59(1):53-8.
6. de Tomasi A. [History of the evolution of diamond burs in dentistry]. [article in Italian] *Odontostomatol Implantoprosesi*. 1976;2(2):72-4.
7. Walsh JP. Critical review of cutting instruments in cavity preparation: I. diamond stone. *Int Dent J*. 1953;4(1):36-43.
8. Siegel SC, Von Fraunhofer JA. Dental cutting: The historical development of diamond burs. *J Am Dent Assoc*. 1998;129(6):740-5.
9. Sharma S, Shankar R, Srinivas K. An epidemiological study on the selection, usage and disposal of dental burs among the dental practioners. *J Clin Diagn Res*. 2014;8(1):250-4.
10. Pilcher ES, Tietge JD, Draughn RA. Comparison of cutting rates among single-patient-use and multiple-patient-use diamond burs. *J Prosthodont*. 2000;9(2):66-70.
11. Siegel SC, von Fraunhofer JA. Assessing the cutting efficiency of dental diamond burs. *J Am Dent Assoc*. 1996;127(6):763-72.
12. Emir F, Ayyildiz S, Sahin C. What is the changing frequency of diamond burs? *J Adv Prosthodont*. 2018;10(2):93-100.
13. Regev M, Judes H, Ben-Hanan U. Wear mechanisms of diamond coated dental burs. *Tribology—Materials Surfaces & Interfaces*. 2010;4(1):38-42.
14. Christensen, GJ. Single-use vs. multi-use diamond rotary instruments. *Clinicians Report*. 2001;25(9):1,3-4.
15. Imran A, Imran H, Ashley MP. Straight to the point: Considering sharp safety in dentistry. *Br Dent J*. 2018;225(5):391-4.
16. Woods R, Amerena V, David P, et al. Sterilisation: Part 1. Instrument preparation. *FDI World*. 1996;5(2):7-10.
17. Morrison A, Conrod S. Dental burs and endodontic files: Are routine sterilization procedures effective? *J Can Dent Assoc*. 2009;75(1):39.
18. Parker MT, Ball LC. Streptococci and aerococci associated with systemic infection in man. *J Med Microbiol*. 1976;9(3):275-302.
19. Whitworth CL, Martin MV, Gallagher M, Worthington HV. A comparison of decontamination methods used for dental burs. *Br Dent J*. 2004;197(10):635-40.
20. Sushma B, Gugwad S, Pavaskar R, et al. Prions in dentistry: A need to be concerned and known. *J Oral Maxillofac Pathol*. 2016;20(1):111-4.
21. Taylor DM. Inactivation of transmissible degenerative encephalopathy agents: A review. *Vet J*. 2000;159(1):10-7.
22. Walker JT, Dickinson J, Sutton JM, et al. Cleanability of dental instruments—Implications of residual protein and risks from Creutzfeldt-Jakob disease. *Br Dent J*. 2007;203(7):395-401.
23. International Organization for Standardization. ISO 11137-3:2017. Sterilization of health care products—Radiation—Part 3: Guidance on dosimetric aspects of development, validation and routine control. June 2017. [iso.org/standard/63841.html](http://iso.org/standard/63841.html)
24. Brinston RM, Wilson BK. Converting to gamma-radiation sterilization: An overview for medical device manufacturers. *Med Device Technol*. 1993;4(4):18-22.
25. Centers for Disease Control and Prevention. Single-use (disposable) devices. April 27, 2018. [cdc.gov/oralhealth/infectioncontrol/faqs/single-use-devices.html](http://cdc.gov/oralhealth/infectioncontrol/faqs/single-use-devices.html)
26. Centers for Disease Control and Prevention. Guideline for isolation precautions: Preventing transmission of infectious agents in healthcare settings (2007). July 2019. [cdc.gov/infectioncontrol/pdf/guidelines/isolation-guidelines-H.pdf](http://cdc.gov/infectioncontrol/pdf/guidelines/isolation-guidelines-H.pdf)



## CE Quiz

To complete this quiz online and immediately download your CE verification document, visit [www.dentallearning.net/SPU4U](http://www.dentallearning.net/SPU4U), then log in to your account (or register to create an account). Upon completion and passing of the exam, you can immediately download your CE verification document. We accept Visa, Mastercard, Discover, and American Express.

1. Archeological evidence from Tuscany, Italy, has shown that "tool-assisted manipulation" to drill into teeth is dated at \_\_\_\_\_.
  - a. 5,000 years ago
  - b. 10,000 years ago
  - c. 13,000 years ago
  - d. 20,000 years ago
2. Hand-held drills could produce how many rotations per minute?
  - a. 15
  - b. 20
  - c. 100
  - d. 60
3. Who invented the first electric drill?
  - a. James Morrison
  - b. George Green
  - c. George Fellows Harrington
  - d. George Harrison
4. Modern methods of removing dental tissue include all of the following except \_\_\_\_\_.
  - a. Lasers
  - b. Air abrasion
  - c. Chemical dissolution
  - d. Radio frequency
5. Which organization oversees the 15-digit code numbering system for dental burs?
  - a. European Union
  - b. International Organization for Standardization
  - c. United Nations
  - d. World Health Organization
6. Burs currently are made from all of the following materials except \_\_\_\_\_.
  - a. Plastic
  - b. Tungsten
  - c. Diamond
  - d. Composite
7. Tungsten carbide is \_\_\_\_ as stiff as steel.
  - a. 1/2 x
  - b. 2 x
  - c. 2.5 x
  - d. 3 x
8. The first diamond bur was invented in \_\_\_\_\_.
  - a. 1930
  - b. 1880
  - c. 1897
  - d. 1902
9. The modern diamond bur was invented in \_\_\_\_\_.
  - a. 1904
  - b. 1924
  - c. 1932
  - d. 1945
10. Which of following manufacturing processes is used to create diamond burs?
  - a. Electro-deposition
  - b. Microbrazing
  - c. Sintering
  - d. All of the above
11. In a survey of dentists, what percentage of respondents said they sterilized burs for every patient?
  - a. 15%
  - b. 35%
  - c. 75%
  - d. 98%
12. Studies have shown similar initial cutting rates for single-patient-use and multiple-patient-use burs.
  - a. True
  - b. False
13. Researchers have recommended that reusable diamond burs should be discarded after how many crown preps?
  - a. 3
  - b. 5
  - c. 7
  - d. 10
14. What is the dominant mechanism of wear of diamond burs?
  - a. Particle wear-out
  - b. Particle pull-out
  - c. Clogging by debris
  - d. Degradation of binder material
15. There is a significant difference in efficiency of single- vs multiple-patient-use diamond burs when judged by time to prepare a lower first molar.
  - a. True
  - b. False
16. In a survey of 164 dentists, what percentage reported having had a sharps injury thus far in their career?
  - a. 5%
  - b. 31%
  - c. 66%
  - d. 89%



17. The most common sharps injuries are due to \_\_\_\_.
- Burs
  - Needles
  - Cavitrons
  - Scalpels
18. The three main methods of sterilization include all of the following, except \_\_\_\_.
- Steam and pressure
  - Dry heat
  - Chemical vapor
  - Enzymatic bath
19. Autoclave alone is sufficient to sterilize burs and endodontic files.
- True
  - False
20. Which organism has been identified as a cause of endocarditis and meningitis?
- Coronavirus
  - Streptococcus sanguinis*
  - Streptococcus mutans*
  - Staphylococcus aureus*
21. According to research, what is the most effective presterilization method to remove debris?
- Manual scrubbing
  - Cavitron bath
  - Enzymatic bath
  - Washer-disinfector
22. Prions contain \_\_\_\_.
- DNA
  - RNA
  - Protein
  - Lipids
23. Which sterilization method is most effective against prions?
- Steam autoclave
  - Enzymatic bath
  - Sodium hypochlorite
  - Dry heat
24. Which country has mandated single-patient-use endodontic files?
- United Kingdom
  - South Korea
  - United States
  - Germany
25. A common manufacturing method used to sterilize burs is \_\_\_\_.
- Enzymatic bath
  - UV-C
  - Dry heat
  - Gamma radiation
26. The Centers for Disease Control and Prevention considers all diamond-coated burs and scaler tips single-patient-use unless what happens?
- Manufacturers submit a 510(k)
  - Manufacturers presterilize each bur
  - Manufacturers can prove the bur won't break
  - Manufacturers pay regulator fees
27. Prions have an incubation period of \_\_\_\_.
- 5 days
  - 1 week
  - 3 months
  - 2-35 years
28. How many people are occupationally exposed to tungsten each year?
- 17,000
  - 27,000
  - 47,000
  - 147,000
29. Diamond burs should only be used once if \_\_\_\_.
- You are prepping a tooth
  - The patient has an infectious disease
  - Doing try-in procedures of metal or zirconia substructure
  - You don't have an autoclave
30. Which of the following is important when considering whether to adopt single-patient-use armamentarium?
- Cost
  - Efficiency
  - Infection control
  - All of the above



**CE ANSWER FORM** (E-mail address required for processing)

WWW.DENTALLEARNING.NET/SPU4U

Name:	Title:	Specialty:
Address:	NPI No.:	
City:	State:	Zip:
AGD Identification No.:		
Email:		
Telephone:	License Renewal Date:	

**AGD Codes: 148, 250**

**EDUCATIONAL OBJECTIVES**

- Discuss the history of dental burs and handpieces
- List the pros and cons of single-patient-use and multiple-patient-use armamentarium
- Analyze dental bur cost, efficiency, and infection control.

**COURSE EVALUATION**

Please evaluate this course using a scale of 1 to 5, where 1 is poor and 5 is excellent.

- Clarity of objectives ..... ① ② ③ ④ ⑤
- Usefulness of content ..... ① ② ③ ④ ⑤
- Benefit to your clinical practice..... ① ② ③ ④ ⑤
- Usefulness of the references ..... ① ② ③ ④ ⑤
- Quality of written presentation ..... ① ② ③ ④ ⑤
- Quality of illustrations ..... ① ② ③ ④ ⑤
- Clarity of quiz questions ..... ① ② ③ ④ ⑤
- Relevance of quiz questions ..... ① ② ③ ④ ⑤
- Rate your overall satisfaction with this course ..... ① ② ③ ④ ⑤
- Did this lesson achieve its educational objectives?  Yes  No
- Are there any other topics you would like to see presented in the future? \_\_\_\_\_
- Overall administration of the program..... ① ② ③ ④ ⑤

**COURSE SUBMISSION:**

- Read the entire course.
- Complete this entire answer sheet in either pen or pencil.
- Mark only one answer for each question.
- Mail or fax answer form to 732-303-0555.

**For immediate results:**

- Read the entire course.
- Go to [www.dentallearning.net/SPU4U](http://www.dentallearning.net/SPU4U).
- Log in to your account or register to create an account.
- Complete course and submit for grading to receive your CE verification certificate.

A minimum score of 70% is required to receive CE credits.

Dental Learning, LLC  
500 Craig Road, First Floor  
Manalapan, NJ 07726

\*If paying by credit card, please note:  
Mastercard | Visa | AmEx | Discover

\*Account Number \_\_\_\_\_

\*Expiration Date \_\_\_\_\_

The charge will appear as *Dental Learning, LLC*.

If paying by check, make check payable to *Dental Learning, LLC*.

ALL FIELDS MARKED WITH AN ASTERISK (\*) ARE REQUIRED

**QUIZ ANSWERS**

Fill in the circle of the appropriate answer that corresponds to the question on previous pages.

- |                         |                         |
|-------------------------|-------------------------|
| 1. (A) (B) (C) (D) (E)  | 16. (A) (B) (C) (D) (E) |
| 2. (A) (B) (C) (D) (E)  | 17. (A) (B) (C) (D) (E) |
| 3. (A) (B) (C) (D) (E)  | 18. (A) (B) (C) (D) (E) |
| 4. (A) (B) (C) (D) (E)  | 19. (A) (B) (C) (D) (E) |
| 5. (A) (B) (C) (D) (E)  | 20. (A) (B) (C) (D) (E) |
| 6. (A) (B) (C) (D) (E)  | 21. (A) (B) (C) (D) (E) |
| 7. (A) (B) (C) (D) (E)  | 22. (A) (B) (C) (D) (E) |
| 8. (A) (B) (C) (D) (E)  | 23. (A) (B) (C) (D) (E) |
| 9. (A) (B) (C) (D) (E)  | 24. (A) (B) (C) (D) (E) |
| 10. (A) (B) (C) (D) (E) | 25. (A) (B) (C) (D) (E) |
| 11. (A) (B) (C) (D) (E) | 26. (A) (B) (C) (D) (E) |
| 12. (A) (B) (C) (D) (E) | 27. (A) (B) (C) (D) (E) |
| 13. (A) (B) (C) (D) (E) | 28. (A) (B) (C) (D) (E) |
| 14. (A) (B) (C) (D) (E) | 29. (A) (B) (C) (D) (E) |
| 15. (A) (B) (C) (D) (E) | 30. (A) (B) (C) (D) (E) |

**Price: \$29 CE Credits: 2**  
Save time and the environment by taking this course online.

**Coupon Code:** \_\_\_\_\_

**If you have any questions, please email [Dental Learning at questions@dentallearning.net](mailto:questions@dentallearning.net) or call 888-724-5230.**

**PLEASE PHOTOCOPY ANSWER SHEET FOR ADDITIONAL PARTICIPANTS.**

Please direct all questions pertaining to Dental Learning, LLC or the administration of this course to [Support@dentallearning.net](mailto:Support@dentallearning.net). COURSE EVALUATION and PARTICIPANT FEEDBACK: We encourage participant feedback pertaining to all courses. Please be sure to complete the evaluation included with the course. INSTRUCTIONS: All questions have only one answer. Participants will receive confirmation of passing by receipt of a verification certificate. Verification certificates will be processed within two weeks after submitting a completed examination. EDUCATIONAL DISCLAIMER: The content in this course is derived from current information and research based evidence. Any opinions of efficacy or perceived value of any products mentioned in this course and expressed herein are those of the author(s) of the course and do not necessarily reflect those of Dental Learning. Completing a single continuing education course does not provide enough information to make the participant an expert in the field related to the course topic. It is a combination of many educational courses and clinical experience that allows the participant to develop skills and expertise. COURSE CREDITS/COST: All participants scoring at least 70% on the examination will receive a CE verification certificate. Dental Learning, LLC is an ADA CERP recognized provider. Dental Learning, LLC is also designated as an Approved PACE Program Provider by the Academy of General Dentistry. The formal continuing education programs of this program provider are accepted by AGD for Fellowship, Mastership, and membership maintenance credit. Please contact Dental Learning, LLC for current terms of acceptance. Participants are urged to contact their state dental boards for continuing education requirements. Dental Learning, LLC is a California Provider. The California Provider number is RP5062. The cost for courses ranges from \$19.00 to \$90.00. RECORD KEEPING: Dental Learning, LLC maintains records of your successful completion of any exam. Please contact our offices for a copy of your continuing education credits report. This report, which will list all credits earned to date, will be generated and mailed to you within five business days of request. Dental Learning, LLC maintains verification records for a minimum of seven years. CANCELLATION/REFUND POLICY: Any participant who is not 100% satisfied with this course can request a full refund by contacting Dental Learning, LLC in writing or by calling 1-888-724-5230. Go Green, Go Online to [www.dentallearning.net](http://www.dentallearning.net) to take this course. © 2020