



Order online: BiteChek.com  
 By phone: 800.235.1863

By fax: 770.423.4996

- Box will dispense small quantities for ease of use.
- Store at room temperature.
- Single-patient-use.



BC-Black  
 Double-Sided  
 19 microns



BC-Red/Black  
 Double-Sided  
 19 microns



## New Dispenser

At Microcopy, we are always striving to make dentistry easier. The new box design dispenses Bite-Chek without risking contamination to the rest of the box during the procedure. This new dispenser makes it easy to protect the box contents and your patient.



	COLOR FILM	100 PER PACK	TOTAL
\$17 <sup>95</sup> per Pack	BC-Black	@ \$17 <sup>95</sup>	\$ _____
	BC-Red/Black	@ \$17 <sup>95</sup>	\$ _____
	US Shipping		\$ <u>6<sup>95</sup></u>
	(Add sales tax if delivered in GA) <b>Total</b>		\$ _____
	Canada Shipping		\$ <u>20<sup>00</sup></u>
	Canada Only - Free Shipping on Orders Over \$400		
	Canadian Orders Apply GST and / or HST		\$ _____
	(Canadian Orders ship direct from Canada) <b>Total</b>		\$ _____



Perry L. Parke  
 President/Owner

### "No Questions Asked Guarantee."

"We unconditionally guarantee your total satisfaction with the quality and performance of all our products. If you are not completely satisfied, we will pick up the product and gladly refund your full purchase price."

Microcopy  
 800.235.1863 | P.O. Box 2017 | Kennesaw, GA 30156-9017

facebook.com/microcopydental @microcopydental

Prices effective November 1, 2016

BC1015

# BITE-CHEK<sup>®</sup>

## Get a Handle on Articulation<sup>®</sup>

"Easy grip" handle eliminates forceps.

Check the bite at every chair with a BITE-CHEK on every tray.



"BITE-CHEK is a wonderful new product which has simplified articulation immensely."

R. King DDS, Grayslake, IL

**NEW**  
 RED/BLACK FILM

800.235.1863  
 BITECHEK.COM



# BITE-CHEK®

## New Color!



Double-Sided  
Red/Black Film  
19 microns

## Family History



Microcopy defined innovation when the first batch of Insta-Neg and Insta-Fix solutions for rapid film processing was bottled back in 1970. Our company continued that winning streak with the introduction of the market-leading NeoDiamond burs. Follow-through products included NeoBurr carbides, NeoDrys super absorbents, Gazelle and Minnow nanocomposite polishers. Quality and patient-safety have been the hallmarks of Microcopy since we began.

There's no question – Microcopy products and your practice are a winning combination.

## Always Ready

For many dental visits the first stop is with the hygienist. During this cleaning there may be sealants, complaints about tooth pain, or an implant may need to be checked for contact pressure. These may require the hygienist or doctor to check the occlusion. Keeping forceps by the chair many times is impractical, so the hygienist must leave to retrieve them, all the while the patient and the doctor are waiting.


Occlusion is fundamental to dental procedures and maintenance, so why not make it more efficient? BITE-CHEK — the perfect assistant for checking occlusion.

- ✓ No mess
- ✓ Cut down on sterilizing forceps and scissors
- ✓ Hygienist Friendly — keep at all stations




Ready to use easy-grip handle requires no forceps

No curling, handle holds film flat



## Adjustable Full Arch

With the patented easy-grip handle, Bite-Chek is also used to mark the full arch by holding two side-by-side as shown. The angle can be adjusted to fit the width of the arch.



## Our Articulation Film Advantage

BITE-CHEK accurately marks the points of contact, minimizing the possibility of “too high” an occlusion and possible post-operative tooth sensitivity. Contrary to many thicker paper options on the market today, our film’s thin profile helps patients close their bite completely, eliminating mandibular reflexes that can skew the occlusion. The result is a small mark that makes adjustment accurate. BITE-CHEK with its easy-grip handle is the ideal tool for precise occlusal adjustments.




Patented  
USD733302 S1


## Bite-Chek Articulating Film:

- ✓ Thin film marks high points accurately – works best on a dry surface.
- ✓ Stretches without snap back or tearing to prevent false marks.
- ✓ Marks all surfaces: enamel, metal, non-metal restorations, and plastic.
- ✓ Is available in Black/Black and Red/Black film.
- ✓ Marks both arches at once with its double-sided film.


## Instructions



Dry the tooth surface.



Position so film is on the teeth facing the center of the mouth.



Check occlusion.