



Order Online:
GazellePolisher.com

By phone:
800.235.1863

By fax:
770.423.4996



Shapes match natural anatomy



Actual Size

| GAZELLE SATIN | | 25 PER PACK | TOTAL |
|--|------------------|---------------|-----------------|
| \$235 PER POLISHER | COMP-CUP 1 | ___ @ \$58.75 | \$ _____ |
| | COMP-DISC 1 | ___ @ \$58.75 | \$ _____ |
| | COMP-POINT 1 | ___ @ \$58.75 | \$ _____ |
| GAZELLE HI-GLOSS | | 10 PER PACK | TOTAL |
| \$445 PER POLISHER | COMP-CUP 2 | ___ @ \$44.50 | \$ _____ |
| | COMP-DISC 2 | ___ @ \$44.50 | \$ _____ |
| | COMP-POINT 2 | ___ @ \$44.50 | \$ _____ |
| NEW! GAZELLE COMBO PACK* | | 10 PER PACK* | TOTAL |
| 7 at \$2.85 SATIN POLISHERS* | COMP-CUP-COMBO | ___ @ \$34.50 | \$ _____ |
| | COMP-DISC-COMBO | ___ @ \$34.50 | \$ _____ |
| 3 at \$4.85 HI-GLOSS POLISHERS* | COMP-POINT-COMBO | ___ @ \$34.50 | \$ _____ |
| *SAME SHAPE PER COMBO PACK | | US Shipping | \$ <u>6.95</u> |
| (Add sales tax if delivered in GA) | | | \$ _____ |
| Total | | | \$ _____ |
| Canada Shipping | | | \$ <u>20.00</u> |
| Canada Only - Free Shipping on Orders Over \$400 | | | \$ _____ |
| Canadian Orders Apply GST and / or HST | | | \$ _____ |
| (Canadian Orders ship direct from Canada) Total | | | \$ _____ |



Go the distance[®]..

NEW! Durable nanocomposite polisher that does not crumble.



“The GAZELLE Hi-Gloss polisher is the best I have seen in 16 years of practice. The gloss is unbelievable.”

S. Grady DDS, Duluth, GA



Perry L. Parke
President/Owner

No Questions Asked Guarantee

“We unconditionally guarantee your total satisfaction with the quality and performance of all our products. If you are not completely satisfied, we will pick up the product and gladly refund your full purchase price.”

NEW COMBO PACK
CONTAINS SATIN & HI-GLOSS
GAZELLEPOLISHER.COM



Microcopy
800.235.1863 | P.O. Box 2017 | Kennesaw, GA 30156-9017

facebook.com/microcopydental @microcopydental

Prices effective January 1, 2015
GZDT 0315



SILITEC PROCESS

Gives the polisher durability to retain its shape and flexibility to refine the tooth surface.

Stop the Crumble

After prepping the tooth and carefully finishing the composite, another successful restoration is almost complete. The front office just notified you that the next patient has arrived, making you keenly aware of the small amount of time left to stay on schedule. You slide the polisher into the handpiece to put the last touch on your work. With the whine of the drill, you see the polisher disintegrate and break into crumbs. So much for a smooth transition to the next patient. This has been the experience for you and dentists everywhere for years. Why does this happen?



Many competitive polishers crumble because nanocomposite has additional glass to increase strength. However GAZELLE is strong enough to polish nanocomposite without falling apart. Like a gazelle, it is a perfect balance between flexibility with durability. Our Silitec™ process infuses the polisher with strength to retain its shape and not crumble while providing the highest gloss value to the restoration.

On your next restoration, choose the polisher that provides the durability, flexibility, gloss, and convenience. Choose GAZELLE and “Go the distance...”

NEW! Combo Pack

Gazelle is now offered in a combo pack with both Satin and Hi-Gloss polishers. This convenient 10-pack contains the same shape in a choice of cups, points or discs. With this assortment of 7 Satin and 3 Hi-Gloss polishers, you can achieve one of the highest gloss values for polishing composites.



GAZELLE is highest rated by these organizations.



Safety and Convenience

GAZELLE is the first and only sterilized polisher. Each polisher is individually packaged and recommended for single-patient-use.



GAZELLE Satin

The Satin polisher is a clinical workhorse for **everyday nanocomposite restoration**. A quick solution to leave a smooth, satin surface, utilizing aluminum oxide in the SILITEC process.

- ◆ Removes **surface striations** and anatomical defects.
- ◆ Leaves a smooth, **satin gloss** surface.
- ◆ **Strong and durable** – will not break down during use.
- ◆ Individually **sterilized** for single-patient-use.



GAZELLE Hi-Gloss


The Hi-Gloss polisher is for the **most demanding nanocomposite esthetic restoration**. A quick solution to leave a smooth, lustrous surface utilizing diamond particles in the Silitec process.

- ◆ Leaves a **liquid luster** gloss.
- ◆ Ideal for **anterior** restorations.
- ◆ **Strong and durable** – will not break down during use.
- ◆ Individually **sterilized** for single-patient-use.




Instructions for Use

- GAZELLE Satin:** apply light pressure.
- GAZELLE Hi-Gloss:** apply feather light pressure.
- Max speed 5,000 rpm.*
- Use no polishing paste.
- Rinse and dry.
- For single-patient-use.



Satin
Light Pressure



Hi-Gloss
Feather Light Pressure

**Avoid long periods of applied pressure as heat generation can cause damage to the margins. Apply water spray to cool if needed.*