

MINNOW®

The same people who make Minnow, also bring you NeoDiamond, NeoBurr, and NeoDrys.



NEODIAMOND®
#1 Single-Patient-Use Diamond



NEOBURR®
New Definition of Strength



NEODRYS®
Ultimate Parotid Saliva Control

PRICING

Contact us directly or one of our dealer partners for pricing.

Order Online: MinnowPolisher.com By phone: 800.235.1863 By fax: 770.423.4996

Minnow Satin

20 pack



MCOMP-CUP 1



MCOMP-POINT 1
(Shapes shown actual size)



MCOMP-COMBO 1

**Combo Packs have 10 Polishers each of the Cup and Point shape.*

Minnow Hi-Gloss

20 pack



MCOMP-CUP 2



MCOMP-POINT 2
(Shapes shown actual size)



MCOMP-COMBO 2

**Combo Packs have 10 Polishers each of the Cup and Point shape.*

MINNOW®

The Go-Anywhere Polisher®

The right fit for the right price.



“ Like the shape and ability to get into grooves and big curves on lingual/buccal surfaces. ”

Dr. Gary Radz

MINNOWPOLISHER.COM

MICROCOPY

Microcopy
800.235.1863 | P.O. Box 2017 | Kennesaw, GA 30156-9017

facebook.com/microcopydental [@microcopydental](https://twitter.com/microcopydental) [@microcopydental](https://instagram.com/microcopydental)



At home, in the tight crevices of rocks and weeds, minnows swim freely where big fish can't reach them. Being nimble and small has many advantages. This is what inspired Microcopy's Minnow polishers.

Perfectly suited for single-patient-use, Minnow provides just the right fit when accessing pits, fissures and tooth anatomy that other, larger polishers just can't reach.

Small Has its Advantages

The compact design of the point and flare of the cup provides easy access to all areas of composite, even next to soft tissue. Their small size provides an efficient and economical answer for finishing off all classes of restorations using today's nanocomposite.

Cost Efficiency is Key

In today's economy, dental professionals need a true single-patient-use, cost-efficient polisher. Minnow comes pre-sterilized — use it and toss it. Also, toss away all the extra costs that go with cleaning and sterilizing multi-use polishers. With no polishing paste needed, Minnow saves you both time and money.

Superior Polishing is the End Result

This simple polishing system gives you a choice between a 1- or 2- step polishing: 1-step for finishing off posterior restorations and 2-step to bring anteriors to a high luster shine.

- ◆ Miniature shape **fits in pits and fissures**.
- ◆ Ideal for **minimally invasive procedures**.
- ◆ Polish **without crumbling** or falling apart.
- ◆ Truly **economical single-patient-use** polisher.
- ◆ Hi-Gloss polishes to a **super high shine**.
- ◆ Individually **sterilized**, ready to use packets.
- ◆ Cup is ideal for **class V restorations**.



Family History



Microcopy pioneered the development of single-patient-use diamonds. Concern for patient safety from cross contamination pointed the way to invent the highest rated bur, NEODIAMOND®, the market leader for single-patient-use. Follow-through products included NEOBURR® carbides, NEODRYS® super absorbents, and BITE-CHEK® articulation products. Quality and patient-safety have been the hallmarks of Microcopy since we began.

There's no question – Microcopy products and your practice are a winning combination.

Miniature shapes for easy access to all restorations; especially those with deep anatomy, interproximal crevices, and areas of close soft tissue.



STEP 1: MINNOW Satin

The Satin polisher is a clinical workhorse for **everyday restorations**. A quick solution to leave a smooth, satin surface, utilizing aluminum oxide in the SILITEC process.

- ◆ **Removes surface striations** and anatomical defects.
- ◆ Leaves a smooth, **satin gloss** surface.



STEP 2: MINNOW Hi-Gloss

The Hi-Gloss polisher is for the **most demanding esthetic restoration**. A quick solution to leave a smooth, lustrous surface utilizing diamond particles in the SILITEC process.

- ◆ Leaves a **liquid luster** gloss.
- ◆ Ideal for **anterior** restorations.



Instructions for Use

MINNOW Satin: apply light pressure.

MINNOW Hi-Gloss: apply feather light pressure.

Max speed 5,000 rpm.*

Use no polishing paste.

Rinse and dry.

For single-patient-use.

**Avoid long periods of applied pressure as heat generation can cause damage to the margins. Apply water spray to cool if needed.*

See all instructions at microcopydental.com/minnow#instructions.



Satin
Light Pressure



Hi-Gloss
Feather Light Pressure