

Minnow is available convenient 20 packs. (Shapes shown actual size)



INTRODUCTORY PRICE	REGULAR PRICE	MINNOW SATIN	20 PER PACK	TOTAL
95¢ PER POLISHER	\$130 PER POLISHER	MCOMP-CUP 1	___ @ \$1895	\$ _____
		MCOMP-POINT 1	___ @ \$1895	\$ _____
		MCOMP-COMBO 1*	___ @ \$1895	\$ _____
INTRODUCTORY PRICE	REGULAR PRICE	MINNOW HI-GLOSS	20 PER PACK	TOTAL
\$175 PER POLISHER	\$200 PER POLISHER	MCOMP-CUP 2	___ @ \$3495	\$ _____
		MCOMP-POINT 2	___ @ \$3495	\$ _____
		MCOMP-COMBO 2*	___ @ \$3495	\$ _____
<small>*COMBO PACKS HAVE 10 POLISHERS EACH OF THE CUP AND POINT SHAPE.</small>				
US Shipping				\$ <u>695</u>
(Add sales tax if delivered in GA)				\$ _____
Total				\$ _____



"No Questions Asked Guarantee."

"We unconditionally guarantee your total satisfaction with the quality and performance of all our products. If you are not completely satisfied, we will pick up the product and gladly refund your full purchase price."

Perry L. Parke
President/Owner



Order Online:
MinnowPolisher.com

By phone:
800.235.1863

By fax:
770.423.4996



Microcopy
800.235.1863 | P.O. Box 2017 | Kennesaw, GA 30156-9017

facebook.com/microcopydental @microcopydental

Prices effective September 1, 2016

MINNOW™

The Go-Anywhere Polisher™

The right fit for the right price.



“ Like the shape and ability to get into grooves and big curves on lingual/buccal surfaces. ”

Dr. Gary Radz

NEW PRODUCT
Introductory Price As Low As
95¢
MINNOWPOLISHER.COM





At home, in the tight crevices of rocks and weeds, minnows swim freely where big fish can't reach them. Being nimble and small has many advantages. This is what inspired Microcopy's new Minnow polishers.

Perfectly suited for single-patient-use, Minnows provide just the right fit when accessing pits, fissures and tooth anatomy that other, larger polishers just can't reach.

Small Has its Advantages

The compact design of the point and flare of the cup lets you access all areas of composite, even next to soft tissue. Their small size provides an efficient and economical answer for finishing off all classes of restorations using today's nanocomposite.

Cost Efficiency is Key

In today's economy, dental professionals need a true single-patient-use, cost-efficient polisher. Minnow comes pre-sterilized — use it and toss it. And, toss away all the extra costs that go with cleaning and sterilizing multi-use polishers. With no polishing paste needed, Minnow saves you both time and money.

Superior Polishing is the End Result

This simple polishing system gives you a choice between a 1- or 2- step polishing: 1-step for finishing off posterior restorations and 2-step to bring anteriors to a high luster shine

- ◆ Miniature shape **fits in pits and fissures**.
- ◆ Ideal for **minimally invasive procedures**.
- ◆ Polish **without crumbling** or falling apart.
- ◆ Truly **economical single-patient-use** polisher.
- ◆ Hi-Gloss polishes to a **super high shine**.
- ◆ Individually **sterilized**, ready to use packets.



Family History



Microcopy defined innovation when the first batch of Insta-Neg and Insta-Fix solutions for rapid film processing was bottled back in 1970. Our company continued that winning streak with the introduction of the market-leading NeoDiamond burs. Follow-through products included NeoBurr carbides, NeoDrys super absorbents, Gazelle and Minnow nanocomposite polishers, and Bite-Check articulation products. Quality and patient-safety have been the hallmarks of Microcopy since we began.

There's no question – Microcopy products and your practice are a winning combination.

Miniature shapes ease access for all restoration; especially those with deep anatomy, interproximal crevices, and areas of close soft tissue.



STEP 1: MINNOW Satin

The Satin polisher is a clinical workhorse for **everyday restoration**. A quick solution to leave a smooth, satin surface, utilizing aluminum oxide in the SILITEC process.

- ◆ **Removes surface striations** and anatomical defects.
- ◆ Leaves a smooth, **satin gloss** surface.



STEP 2: MINNOW Hi-Gloss (Optional)

The Hi-Gloss polisher is for the **most demanding esthetic restoration**. A quick solution to leave a smooth, lustrous surface utilizing diamond particles in the Silitec process.

- ◆ Leaves a **liquid luster** gloss.
- ◆ Ideal for **anterior** restorations.



Instructions for Use

MINNOW Satin: apply light pressure.

MINNOW Hi-Gloss: apply feather light pressure.

Max speed 5,000 rpm.*

Use no polishing paste.

Rinse and dry.

For single-patient-use.

**Avoid long periods of applied pressure as heat generation can cause damage to the margins. Apply water spray to cool if needed.*



Satin
Light Pressure



Hi-Gloss
Feather Light Pressure